

- Minimum Hardware Requirements:
 - Minimum 4GB Ram
 - 1GHz CPU
 - Minimum i3 processor
 - 1366 x 768 Display Resolution
 - 1 USB Port
 - 1 GB free space
- Minimum Software Requirements:
 - Microsoft Windows8, Windows10, Windows11
 - .Net Framework 4.5 Full

The latest Fitting version below link:

http://dosya.eartechnic.com/EarTuning2/EarTuning2_Setup.rar

http://dosya.eartechnic.com/CorrectFit2/CorrectFit2_Setup.rar

Note1: Unzip the software file before setting up

Note2: Make sure the antivirus and firewall software are turned off during the download of EarTuning2 or CorrectFit2 Software.

Note3: Do not turn off or restart your computer during the download.

Note4:

- HiPro and HiPro2:
 - ✓ Make sure the original HiPro device and software used.
 - ✓ Configure HiPro

PREPARED BY

CONTROLLED BY

APPROVED BY

Production Manager

Quality Assurance Specialist

Quality Management Representative

Requirements for Downloading Fitting Software to Computer

Document No	TD.01.DE.112
Release Date	20.03.2025
Revision Date	13.08.2025
Revision No	02
Page No:	2 / 2

User Roles and Permissions:

The user of the Fitting Software is always the clinical expert.

Clinical experts can access all parameters of the device through the fitting program and make the appropriate adaptations to the patient.

Automatic Software Updates:

Fitting Software must be maintained with continuous security patches and updates.

So Automatic updates or manual updates of the software should be provided.

After the program update is completed, a “What's New” window automatically pops up to inform the user.

This screen provides information about all the new features and changes in bullet points.

PREPARED BY

CONTROLLED BY

APPROVED BY

Production Manager

Quality Assurance Specialist

Quality Management Representative